

## ULTRASONIC DIAGNOSTIC APPARATUS

## ABSTRACT OF THE DISCLOSURE

For the purpose of allowing confirmation of time progress during a run of a  
5 contrast agent protocol and the like on a display screen, time elapsed after  
freeze-off and time duration of freeze are displayed in a contrast agent protocol  
display region 52. In an elapsed-time-after-injection display region 55, time  
elapsed from the time of contrast agent injection is displayed. In a countdown  
display region 56, time remaining up to a scheduled time for the next operation  
10 to be made by an operator is calculated and displayed.

BEST AVAILABLE COPY