

ULTRASONIC DIAGNOSTIC APPARATUS

ABSTRACT OF THE DISCLOSURE

For the purpose of allowing confirmation of time progress during a run of a
5 contrast agent protocol and the like on a display screen, time elapsed after
freeze-off and time duration of freeze are displayed in a contrast agent protocol
display region 52. In an elapsed-time-after-injection display region 55, time
elapsed from the time of contrast agent injection is displayed. In a countdown
display region 56, time remaining up to a scheduled time for the next operation
10 to be made by an operator is calculated and displayed.

BEST AVAILABLE COPY